


Information och måttabell

Rund- & flatstickat, arm standardstorlekar



mediven esprit

Armstrumpor omkrets (cm)

Storlek	I	II	III	IV	V
cA	17-19	18-20	19-22	21-24	23-27
cC	15-17	17-19	19-21	21-23	23-25
cE	25-27,5	27-29,5	29-32	31-34	33-37
cG	28-32	30-34	32-37	35-40	38-44

Standardlängd cG: 40-44cm, lång: 44-48cm

Handskar omkrets (cm)

Storlek	I	II	III	IV	V
A	17,5-19,5	18,5-20	20-22	22-24	24-27
B	18-19	19-21,5	21,5-23,5	23,5-25,5	25,5-28,5
C	15-17	17-19	19-21,5	21-23	23-25

Handskar längd (cm)

AB	3,5-4,5
AC	10-12
AC ₁	15-16

Fingrar längd (cm)

Tumme	3,5-4,5
Pekfinger	4,5-5,5
Mittfinger	5,5-6,5
Ringfinger	4,5-5,5
Lillfinger	3,5-4,5

mediven

Armstrumpor omkrets (cm)

Storlek	II	III	IV	V	VI	VII
cA	16-19	18-21	20-23	22-25	24-27	26-29
cC	13-15	15-17	17-19	19-21	21-23	23-25
cE	20-26	23-29	26-32	29-35	32-38	35-41
cG	22-29	25-32	29-36	33-40	37-44	41-48

Standardlängd cG: 40-44cm, lång: 44-48cm

mediven 95

Armstrumpor omkrets (cm)

Storlek	II	III	IV	V	VI	VII
cA	16-19	18-21	20-23	22-25	24-27	26-29
cC	13-15	15-17	17-19	19-21	21-23	23-25
cE	20-26	23-29	26-32	29-35	32-38	35-41
cG	22-29	25-32	29-36	33-40	37-44	41-48

Standardlängd cG: 40-44cm, lång: 44-48cm

Handskar omkrets (cm)

Storlek	II	III	IV	V	VI	VII
A	17,5-19,5	18,5-20	20-22	22-24	24-27	26-29
B	18-19	19-21,5	21,5-23,5	23,5-25,5	25,5-28,5	26,5-30
C	15-17	17-19	19-21	21-23	23-25	23-25
C ₁	14-16,5	16-18,5	18-20,5	20-22,5	22-24,5	24-26,5

Handskar längd (cm)

AB	3,5-4,5
AC	10-12
AC ₁	15-16

Fingrar längd (cm)

Tumme	3,5-4,5
Pekfinger	4,5-5,5
Mittfinger	5,5-6,5
Ringfinger	4,5-5,5
Lillfinger	3,5-4,5

Fingrar omkrets (cm)

Storlek	II	III	IV	V	VI	VII	
Tumme	Topp	5,5	6,5	7,5	8,5	9,5	10,5
	Bas	6,5	7,5	8,5	9,5	10,5	11,5
Pekfinger	Topp	4,5	5,5	6,5	7,5	8,5	9,5
	Bas	5,5	6,5	7,5	8,5	9,5	10,5
Mittfinger	Topp	4,5	5,5	6,5	7,5	8,5	9,5
	Bas	5,5	6,5	7,5	8,5	9,5	10,5
Ringfinger	Topp	4,0	5,0	6,0	7,0	8,0	9,0
	Bas	5,0	6,0	7,0	8,0	9,0	10,0
Lillfinger	Topp	3,5	4,5	5,5	6,5	7,5	8,5
	Bas	4,5	5,5	6,5	7,5	8,5	9,5