


# Information och måttabell

## Rund- & flatstickat, arm standardstorlekar



### mediven harmony

Armstrumpor omkrets (cm)								Längd (cm)	
Storlek	I	II	III	IV	V	VI	VII	Standardlängd	41 – 47
cG x-vid	28,5-31,0	30,5-33,5	33,0-36,0	35,0-38,5	37,0-40,5	39,5-43,0	41,5-45,5		
cG	25,5-28,0	28,0-30,5	30,0-33,0	32,0-35,0	33,5-36,5	36,0-39,5	38,0-41,5		
cF x-vid	25,5-28,0	27,5-30,0	30,0-33,0	32,0-35,0	34,0-37,0	36,0-39,5	38,5-42,0		
cF	23,5-25,5	25,0-27,5	27,5-30,0	29,5-32,0	31,0-34,0	33,0-36,0	35,0-38,5		
cE x-vid	24,0-26,0	26,5-29,0	28,5-31,0	30,0-33,0	31,5-34,5	33,5-36,5	35,5-39,0		
cE	23,0-25,0	24,5-26,5	26,0-28,5	27,5-30,0	29,0-31,5	30,5-33,5	32,5-35,5		
cD x-vid	24,0-26,0	25,0-27,5	26,0-28,5	28,0-30,5	29,5-32,0	31,0-34,0	33,5-36,5		
cD	22,0-24,0	23,0-25,0	24,0-26,0	25,0-27,5	27,0-29,5	28,5-31,0	30,5-33,5		
cC <sub>1</sub> = 100%*	15,0-16,5	16,5-18,0	18,0-19,5	19,5-21,0	21,0-23,0	23,0-25,0	25,0-27,5		
cC	13,5-14,5	14,5-15,5	16,0-17,5	17,5-19,0	19,0-20,5	20,5-22,5	22,5-24,5		
cA	15,0-18,0	16,0-19,0	17,0-20,0	19,0-22,0	20,0-23,0	22,0-25,0	24,0-27,0		

\* Mätpunkt cC<sub>1</sub> är 5 cm ovanför cC

Handskar omkrets (cm)							Fingrar omkrets (cm)							
Storlek	II	III	IV	V	VI	VII	Storlek	II	III	IV	V	VI	VII	
cC <sub>1</sub> = 100%*	15,0-18,0	17-20	19-21,5	21-23,5	23-26	25,5-30	Tumme	Bas	7,0	7,5	8,0	8,5	9,0	10,5
cC	13,5-15,5	15-18	17-20	19-21	20-23,5	22,5-27,5	Tumme	Topp	6,0	6,5	6,5	7,0	7,5	8,0
cB	15,5-20	17,5-21	19,5-23	20,5-24	22,5-26	24,5-30	Pekfinger	Bas	6,5	7,0	7,5	8,0	8,5	9,5
cA	15-19	17-20	19-20	20-23	22-25	24-29	Pekfinger	Topp	5,0	5,5	6,0	6,0	6,5	7,0
							Mittfinger	Bas	6,5	7,0	7,5	8,0	8,5	9,5
							Mittfinger	Topp	5,0	5,5	6,0	6,0	6,5	7,0
							Ringfinger	Bas	6,0	6,5	7,0	7,0	7,5	9,0
							Ringfinger	Topp	5,0	5,0	5,5	5,5	6,0	6,5
							Lillfinger	Bas	5,5	6,0	6,0	6,5	7,0	8,0
							Lillfinger	Topp	4,5	4,5	5,0	5,5	5,5	6,0

\* Mätpunkt cC<sub>1</sub> är 5 cm ovanför cC

### mediven esprit

Armstrumpor omkrets (cm)						Handskar omkrets (cm)					
Storlek	I	II	III	IV	V	Storlek	I	II	III	IV	V
cA	17-19	18-20	19-22	21-24	23-27	A	17,5-19,5	18,5-20	20-22	22-24	24-27
cC	15-17	17-19	19-21	21-23	23-25	B	18-19	19-21,5	21,5-23,5	23,5-25,5	25,5-28,5
cE	25-27,5	27-29,5	29-32	31-34	33-37	C	15-17	17-19	19-21,5	21-23	23-25
cG	28-32	30-34	32-37	35-40	38-44						

Standardlängd cG: 40-44cm, lång: 44-48cm

Handskar längd (cm)		Fingrar längd (cm)	
AB	3,5-4,5	Tumme	3,5-4,5
AC	10-12	Pekfinger	4,5-5,5
AC <sub>1</sub>	15-16	Mittfinger	5,5-6,5
		Ringfinger	4,5-5,5
		Lillfinger	3,5-4,5