



circaid® juxtafit® Måttabeller för ben

circaid® juxtafit® essentials & premium underben



Omfångsmått angivna i cm	S	M	M - X	L
eC	30 - 38	36 - 43	43 - 50	41 - 48
eB1	20 - 28	25 - 33	33 - 40	30 - 38
eB	15 - 23	20 - 28	20 - 28	25 - 33
	L - X	XL	XL - X	XXL
eC	48 - 56	41 - 51	51 - 61	53 - 64
eB1	38 - 46	33 - 41	41 - 51	43 - 53
eB	25 - 33	30 - 38	30 - 38	33 - 41
Produktens längd i cm		kort	lång	
lB - D		28	36	
Rekommenderat längdmått i cm		kort	lång	
l fotknöl - knäveck		< 36	> 36	


circaid® juxtafit® essentials övre ben med knä



Omfångsmått angivna i cm	XXS	XS	S	M	L	XL
eG	40 - 52	46 - 62	54 - 70	60 - 80	69 - 90	77 - 101
eE1*	32 - 42	37 - 48	42 - 55	50 - 65	58 - 76	66 - 87
eD	23 - 30	27 - 35	31 - 41	35 - 46	40 - 52	44 - 58
Produktens längd i cm				Extra kort	Kort	Lång
lD - ljumske			35	45	55	
Rekommenderat längdmått i cm				Extra kort	Kort	Lång
lE - ljumske			< 25	25 - 35	> 35	

* 5 cm över patella

circaid® fot och vrist



circaid juxtafit premium interlocking afw			circaid power added compression band (pac band)	
	S	M	Standard	Large
lA - Y	> 9	> 9		
eB	< 36	< 36	19 - 27	25 - 50
eH	19 - 24	24 - 29	19 - 36	22 - 41
lA	13 - 25	13 - 25		
circaid juxtafit premium afw			circaid single band afw	
	S	M	Bredd	Längd
lA - Y	> 9	> 9	7,6 cm	203 cm
eB	< 36	< 36		
eH	19 - 24	24 - 29		
lA	> 16	> 16		
			circaid ankelsocka KKL 1	
			Standard	Large
eB			19 - 27	25 - 50
eH			19 - 36	22 - 41